

垃圾袋 Garbage Bags

在医疗机构和实验室，不同颜色和材质的垃圾袋用于处理和回收不同的医疗废弃物。
Different colors and material of the garbage bags are used for different purposes in hospital and laboratory.

HDPE T Shirt Garbage Bags (背心手提式)

Item #	Size(W+G)xL	Trash Can Perimeter	Trash Can Height	Capacity	Case Pack
GB-T-5L	(21+11)x38cm	Sharps Box ≤ 5L(5L 以下利器盒)			50pcs
GB-T-10L	(28+14)x48cm	P ≤ 84cm	H ≤ 23cm	10L	50pcs
GB-T-15L	(35+15)x56cm	P ≤ 100cm	H ≤ 27cm	15L	50pcs
GB-T-30L	(40+20)x70cm	P ≤ 120cm	H ≤ 37cm	30L	50pcs
GB-T-50L	(47+23)x80cm	P ≤ 140cm	H ≤ 43cm	50L	50pcs
GB-T-70L	(55+25)x90cm	P ≤ 160cm	H ≤ 51cm	70L	50pcs

HDPE/LDPE/PP or Biodegradable Material Flat Opening Garbage Bags (平口式)

Item #	Size(W+G)xL	Trash Can Perimeter	Trash Can Height	Capacity	Case Pack
GB-F-20L	50x60cm	P ≤ 100cm	H ≤ 38cm	20L	50pcs
GB-F-30L	60x70cm	P ≤ 120cm	H ≤ 46cm	30L	50pcs
GB-F-50L	70x80cm	P ≤ 140cm	H ≤ 53cm	50L	50pcs
GB-F-70L	80x90cm	P ≤ 160cm	H ≤ 60cm	70L	50pcs
GB-F-90L	80x100cm	P ≤ 160cm	H ≤ 65cm	90L	50pcs
GB-F-100L	90x110cm	P ≤ 180cm	H ≤ 70cm	100L	50pcs
GB-F-120L	(80+20)x120cm	P ≤ 200cm	H ≤ 81cm	120L	50pcs
GB-F-240L	(90+30)x140cm	P ≤ 240cm	H ≤ 95cm	240L	50pcs

HDPE/LDPE or Biodegradable Material Drawstring Garbage Bags (抽绳式)

Item #	Size(W+G)xL	Trash Can Perimeter	Trash Can Height	Capacity	Case Pack
GB-D-20L	52x52cm	P ≤ 104cm	H ≤ 30cm	20L	50pcs
GB-D-30L	65x65cm	P ≤ 130cm	H ≤ 40cm	30L	50pcs
GB-D-50L	73x80cm	P ≤ 146cm	H ≤ 50cm	50L	50pcs
GB-D-70L	82x85cm	P ≤ 164cm	H ≤ 55cm	70L	50pcs
GB-D-90L	82x95cm	P ≤ 164cm	H ≤ 65cm	90L	50pcs
GB-D-120L	100x110cm	P ≤ 200cm	H ≤ 75cm	120L	50pcs

HDPE/LDPE or Biodegradable Material Normal Black Garbage Bags (常规黑色垃圾袋)

Item #	Size(W+G)xL	Trash Can Perimeter	Trash Can Height	Capacity	Case Pack
GB-B-12L	45x50cm	P ≤ 104cm	H ≤ 40cm	12L	50pcs
GB-B-20L	50x60cm	P ≤ 120cm	H ≤ 40cm	20L	50pcs
GB-B-30L	60x70cm	P ≤ 120cm	H ≤ 50cm	30L	50pcs
GB-B-50L	70x80cm	P ≤ 140cm	H ≤ 55cm	50L	50pcs
GB-B-70L	80x90cm	P ≤ 160cm	H ≤ 60cm	70L	50pcs
GB-B-90L	80x100cm	P ≤ 160cm	H ≤ 70cm	90L	50pcs
GB-B-100L	90x110cm	P ≤ 180cm	H ≤ 75cm	100L	50pcs
GB-B-120L	100x120cm	P ≤ 200cm	H ≤ 90cm	120L	50pcs
GB-B-240L	120x140cm	P ≤ 240cm	H ≤ 110cm	240L	50pcs



高温灭菌医疗废物袋 Autoclavable Biohazard Bags

用于感染性医疗废弃物的高温高压灭菌处理，温度达到 275° F (135° C) 时，白色指示标将变为灰色或者黑色，标志着印有生物危害标志袋子里的内容物成功灭菌。

Infectious waste is sterilized after collected into autoclave biohazard bags. The bags are ideal for processing hazardous materials up to 275° F (135° C). Imprinted with a biohazard symbol and an indicator which changes to "Autoclaved" or changes to gray when contents are successfully sterilized.

Item #	Size W x L	Case Pack
AC12X242M	12"x24"	50-200pcs
AC14X192M	14"x19"	50-200pcs
AC19X242M	19"x24"	50-200pcs
AC24X302M	24"x30"	50-200pcs
AC23X352M	25"x35"	50-200pcs
AC31X382M	31"x38"	50-200pcs
AC38X472M	38"x47"	50-200pcs
Custom size, logo, color, thickness and packing available		

耐高温医疗废物袋 Biohazard Bags Pass High temperature resistance 121°C/135°C



医用 CT 袋、X 光胶片袋 X-Ray Film Bag

Item #	X-Ray Film Size	Bag Size	Case Pack
CTB-810	8"x10"	25x35cm	50pcs
CTB-8-A4	A4	25x40cm	50pcs
CTB-8-1012	10"x12"	30x40cm	50pcs
CTB-8-1114	11"x14"	30x45cm	50pcs
CTB-8-1215	12"x15"	35x48cm	50pcs
CTB-8-1417	14"x17"	39x52cm	50pcs
Custom size, logo, color, thickness and packing available			



病服袋 Patient Belongings Bags

用于运送收集医疗中心或医院的病服、床单，常规有透明、白色或蓝色的，款式有抽绳的，硬提手的和贴胶孔等，可以定制印刷内容。
Used for transporting or organizing supplies for healthcare providers or for hospital. Available in color of clear, white or blue imprinted bags with various closures to keep contents secure.

Item #	Size W x L	Style/Color	Case Pack
PBB20202D5	20"x20"	drawstring, white, printed	250pcs
PBB18212D5	18"x21"	drawstring, clear, printed	250pcs
PBB20202D5B	20"x20"	drawstring, blue, printed	250pcs
PBB18212DT	18"x21"	drawtape, white, printed	250pcs
PBB120202DT	20"x20"	drawtape, clear, printed	250pcs
PBB18212RH	18"x21"	rigid handle, white, printed	250pcs
PBB2020TE	20"x20"	patch handle, temper evident, printed	250pcs
Custom size, logo, color, thickness and packing available			





药片袋 / The Pill Bags/Medicine Bags

Item #	Size W x L	Case Pack
PBZZ203015	3"x4"	100pcs
PBZZ20909	9x9cm	100pcs
PBZZ201012	10x12cm	100pcs
Custom size, logo, color, thickness and packing available		



污物袋 / Soiled Linen Bags

Regular and drawstring versions. Printed and non-printed both available.

Item #	Size W x L	Case Pack
SLB293802	29" x38"	250pcs
SLB304302	30"x43"	250pcs
Custom size, logo, color and thickness available		



95kPa Biohazard Specimen Bags / 95kPa 样本运输袋

- Passes transportation requirements for International Air Transport Association (IATA) for shipping biohazardous specimens via air
- Engineered to withstand pressure difference of 95kPa
- Adhesive closure keeps contents secure
- Easy-open bag has 2 tear points

Item #	Size W x L	Description	Case Pack
95KPASB0708315	7"x8"	3-wall, PE/PET, 3.15mil	500-1000pcs
95KPASB1215315	12"x15"	3-wall, PE/PET, PA/PE, 3.15mil	500-1000pcs
Custom size, logo, color, thickness and packing available			



药房打包袋 / Pharmacy Bag

Item #	Size W x L	Case Pack
Paper Bag	6" x 3.6" x 11"	100pcs
Pe Bag-T	(32+14)x50cm	100pcs
Non Woven Bag	23cmx30cm	100pcs
Custom size, thickness, color and packing available		



呕吐袋 / 呕吐包 / Vomit Bag/Emesis Bag

一次性的呕吐袋的防逆流设计，可以有效防漏、防异味
Single use emesis bags are designed to reduce leakage and odor.

Item #	Size W x L	Case Pack
VB-PC-1500ML	15.5x36cm	50pcs
VB-RR-1000ML	16.5x40cm	25pcs
VB-DB-1500ML	25x25cm	20pcs
Custom size, logo, color and thickness available		



Yellow Bag with Adhesive Tape / 黄色贴胶袋

常用于医用手术包里面，用来装手术中的一些污物，一般是单只折叠入手术包。
Used to pack into a surgical kit, to load the operating room dirt.
Usually packed piece by piece.

Item #	Size W x L	Description
YAB11516514ADH	11.5"x16.5"	HDPE bag with adhesive tape, 1.4mil
Custom size, logo, color, thickness and packing available		

琥珀色自封袋 Amber Bags

2mil 厚的密封袋，可以有效防止紫外线对药品的损害。
Protect light sensitive contents from damaging ultraviolet light.

Item #	Size W x L
ABZ20302	2"x3"
ABZ250902	2.5"x9"
ABZ30502	3"x5"
ABZ40602	4"x6"
ABZ801402	8"x14"
ABZ901202	9"x12"
ABZ1201202	12"x12"
Custom size, logo, color, thickness and packing available	



样本袋 Specimen Transport Bags

Specimen bag improves lab processing efficiency and accuracy - saving time and improving care.
Dual pouch design features a reclosable poly bag for secondary leakage containment plus an attached open-end document pouch to keep paperwork with the specimen.

Item #	Size W x L	Description	Case Pack
STB60902	6"x9"	4-wall, red & black	1000pcs
STB601002	6"x10"	3-wall, orange & black	1000pcs
STB60102ADH	6"x10"	3-wall, red & black, adhesive tape seal	1000pcs
STB80102	8"x10"	3-wall, orange & black	1000pcs
STB12152	12"x15"	3-wall, orange & black	1000pcs
STB1215R	12"x15"	2-wall, red opaque (no pocket)	1000pcs
Custom size, logo, color, thickness and packing available			



水溶洗衣袋 PVA Water-Soluble Laundry Bag

水溶性 PVA 薄膜是一种新型绿色环保材料，它是一种无污染、无毒的水溶性薄膜，遇水即溶，其水溶液渗入土壤中后 PVOH 高分子易被菌株分解，最终降解为水和二氧化碳 (CO₂ + H₂O)。在医疗机构常用特定的支架收集好医院的病人用过的服装及床单、被罩等脏衣物。半水溶、全水溶材质均可定制。

LOW TEMPERATURE (≥ 25° C)
MIDDLE TEMPERATURE (≥ 45° C)
HIGH TEMPERATURE (≥ 65° C)

Can be decomposed into water and carbon dioxide under the action of natural environment, no harm to the environment.
PVA anti-infection laundry bag has the property of water solution.
After the hospital staff put dirty clothes into the bag, they can be directly put into the washing machine for high-temperature washing and disinfection. The bag will dissolve completely in water during washing.

Item #	Size W x L	Case Pack
WSB26331M	26"x33"	250pcs
WSB28391M	28"x39"	250pcs
WSB36391M	36"x39"	250pcs
Custom size, logo, color, thickness and packing available		



透明和白色穿绳袋 Clear and White Drawstring Bags

两边抽绳的设计可以有效保护袋子内的物品。
Double drawstring closure secures content.

Item #	Size W x L	Case Pack
CDB901215M	9"x12"	250pcs
CDB121815M	12"x18"	250pcs
CDB2024-4M	20"x24"x4"	250pcs
Custom size, logo, color, thickness and packing available		

坐便器内衬袋 (马桶清洁袋) Bedside Commode Liner

The commode liners are universally sized to fit both standard and bariatric commode pails.
Reducing splashes and spills, each leak-proof commode liner includes a super absorbent pad that traps unpleasant odors and holds up to 16 ounces.
An integrated drawstring also allows for quick, spill-proof disposal.

Item #	Size W x L	Case Pack
CL142001	14.76"x20.67"	20pcs/box
Custom size, logo, color, thickness and packing available		



仪器罩 Equipment Covers

一次性的罩子使医疗设备保持清洁，颜色可以定制，不同颜色的罩子便于区分清洁或污染的设备，口部可加皮筋。
Durable pe covers protect and keep medical equipment clean.

一次性床罩 Bed Cover

常见的一次性床罩有 PE/CPE 的，SMS 或者淋膜无纺布，颜色可以定制，口部可加皮筋。
Use material like PE/ CPE/ SMS/ laminated non woven, blue or any other color as your requirement, normally used in hospital, beauty salon, etc. With elastic.

Item #	Size W x L	Case Pack
BC102190	102cmx190cm	100pcs
BC140240	140cmx240cm	100pcs
Custom size, logo, color, thickness and packing available		

枕头罩 Pillow Cover

无纺布枕头罩，颜色可以定制，常用于酒店，医院，美容店等，可以做加皮筋口的，也可以做敞口的，包装方式可定制。
Non-woven material, white, blue or any color as your requirement. Normally used in hotel, hospital, beauty salon, household, etc.
With elastic opening/ square ends or plain.

Item #	Size W x L	Case Pack
PC6060	60cmx60cm+10cm	100pcs
Custom size, thickness, color and packing available		

急救毯 First Aid Blanket

- PET 材质的急救毯，常用于户外急救，可以防风、防水、保暖，常见有银色和金色两种。
- Also called emergency blanket.
- Material: PET
- Color: Gold Silver/Silver Silver
- Feature: Wind proof, waterproof & against coldness



Item #	Size W x L	Case Pack
SSCB3610	36"x10"	250pcs
SSCB2556	25.5"x6"	250pcs
Custom size, logo, color, thickness and packing available		

警示标签 Warning Labels

亮色的医疗警示标签，可以有效的提高使用者的安全警惕。

Bright colored labels improve patient safety. Our labels can be made both in paper and plastic, both on roll and piece by piece are ok.



Item #	Model	Size W x L	Labels/Roll
NS33	No Smoking	3"x3"	100pcs
PS21	Pressure Set	2"x1"	100pcs
CD21	Concentrator Data	2.5"x1.5"	100pcs
REF43	Refrigerate	4.75"x3"	100pcs
DND33	Do Not Discard	3"x3"	100pcs
Custom size, logo, color, thickness and packing available			

急救标识片材 First Aid Plastic Sheet

LDPE 塑料急救片材常用于户外自驾游，飞行，划船，徒步旅行出现意外时。

LDPE plastic emergency sheets are often used for outdoor driving trips, flying, boating, hiking accidents...



Item #	Size W x L	Case Pack
HELP	36"x15"	1pc/bag
SOS	36"x15"	1pc/bag
Custom size, thickness and packing available		

医用计数袋 Surgical Sponge Counter Bag

医用计数袋可以有效解决医用纱布清点不便的问题，袋子有 10 个清晰的口袋，袋子顶部设计有两个挂孔，可以悬挂在合适的架子上。

Surgical Sponge Counter Bag is designed to streamline the process of tracking used surgical sponges in operating room settings by providing a convenient sponge receptacle that also aids in counting.



Item #	Size W x L	Case Pack
SSCB3610	36"x10"	250pcs
SSCB2556	25.5"x6"	250pcs
Custom size, logo, color, thickness and packing available		

信步全自动鞋套机 (Non Power) CHNPACK Automatic Shoe Cover Dispenser

CS100 鞋套机一次可装鞋套 100 只, CS30 鞋套机一次可装鞋套 30 只。使用寿命 10 年以上, 每台鞋套机出厂前都会配一把鞋套。

CS100 can load 100pcs shoe covers per time; CS30 can load 30pcs shoe covers per time. 100% Correct Wearing Rate, 10 years working life. Each dispenser has been installed full cassette PE shoe covers inside.

Item #	Material	Weight	Size W x L	Color	Mode	Case Pack
CS100	ABS+PVC	5.5KGS	23x65x95cm	Silver	non power	1set
CS30	ABS+PVC	2.9KGS	58x23x18cm	Silver, Blue, Red, Yellow, Green	non power	4sets

信步机用鞋套 Shoe Covers For Dispenser

Item #	Weight	Size W x L	Case Pack
CS100HDPE	0.01mm	15cmx50cm	200pcs
CS100CPE	4g/pc	15cmx50cm	200pcs
CS100N-W	23gsm	15cmx50cm	200pcs
CS30HDPE	0.01mm	15cmx50cm	60pcs
CS30CPE	4g/pc	15cmx50cm	60pcs

一次性鞋套 / 靴套 Shoe Covers and Boot Covers

一次性鞋套可以做 HDPE/LDPE/CPE/ 无纺布 / 淋膜无纺布等材质
LDPE/ HDPE/CPE and Non Woven are available.

Item #	Size W x L	Case Pack
SC153B	15cmx38cm	100pcs
SC1540	15cmx40cm	100pcs
SC1550	15cmx50cm	100pcs

Custom size, logo, color, thickness and packing available

无纺布帽 Non Woven Cap

Item #	Thickness	Size W x L	Case Pack
NWC18	10gsm	18"	100pcs
NWC19	12gsm	19"	100pcs
NWC21	15gsm	21"	100pcs

Custom size, logo, color, thickness and packing available



应急尿袋 Emergency Urine Bag



Item #	Size W x L	Case Pack
EUB600CC	13x26.5cm	100pcs
EUB700CC	13x31cm	100pcs
EUB800CC	13x34cm	100pcs

Custom size, thickness, color and packing available

一次性袖套 PE Sleeve Cover



Item #	Size W x L	Case Pack
PSC2040	20cmx40cm	100pcs
PSC2046	20cmx46cm	100pcs

Custom size, thickness, color and packing available

一次性手套 Pe Gloves



Item #	Size W x L	Case Pack
S	27cmx24cm	100pcs
M	28cmx26cm	100pcs
L	29cmx26cm	100pcs

Custom size, thickness, color and packing available

一次性围裙 PE Apron



Item #	Size W x L	Case Pack
PA2442	24"x42"	100pcs
PA2846	28"x46"	100pcs
PA31549	31.5"x49"	100pcs

Custom size, thickness, color and packing available

一次性半椅套/全椅套保护膜 Disposable Half/Full Chair Sleeves

Item #	Size W x L	Case Pack
UD-B11	Full Chair (29"x80")	125pcs/box
UD-B11SG	S&G Full Chair (29"x80")	125pcs/box
UD-B12	Full Chair (29"x80")	125pcs/box

Custom size, thickness, color and packing available



一次性牙椅头靠保护套 Disposable Headrest Cover Sleeves

Item #	Size W x L	Case Pack
UD-B13	Small (11"x9 1/2"x2")	250pcs/box
UD-B13-B	Small (11"x9 1/2"x2")	12x250pcs/box
UD-B14	Large (14"x9 1/2"x2")	250pcs/box
UD-B14-B	Large (14"x9 1/2"x2")	12x250pcs/box

Custom size, thickness, color and packing available



一次性托盘保护膜 Disposable Tray Sleeves

Item #	Size W x L	Case Pack
UD-B01-E	For E-Tray (11 1/2"x16")	500pcs/box
UD-B01-A	For A-Tray (11 5/8"x14 1/2")	500pcs/box
UD-B01-B	For B-Tray (10 1/2"x14")	500pcs/box
UD-B01-F	For F-Tray (7 1/2"x10 1/2")	500pcs/box

Custom size, thickness, color and packing available



一次性键盘保护套 Disposable Keyboard Sleeves

Item #	Size W x L	Case Pack
UD-B58	Small (18.5x55cm)	500pcs/box
UD-B08	Medium (20.5x65cm)	250pcs/box
UD-B57	Large (35.5x56cm)	250pcs/box

Custom size, thickness, color and packing available

